

Face Shields

Enclosed are face shields from the Maker Nexus community of makers. These shields are provided free of charge as humanitarian aid.

We are a network of community members who make things. We are NOT a medical manufacturing company. As such, you are responsible for ascertaining that these face shields are suitable for your use. Maker Nexus and our extended community of makers offer no warranty and accept no liability whatsoever in relation to your use of these masks. If you are unwilling to accept complete and full liability for your use of these face shields, then you should either return them to us or dispose of them as you see fit. We are sorry to have to be so blunt, but this is a reality of the world we live in.

This shipment is NOT sterile. You must take whatever steps you deem necessary to prepare them before use.

Materials

The shield part is made of PETG and is resistant to alcohol and UV light.

The head band and stabilizer could be made of PETG, ABS, or PLA. All materials are resistant to alcohol. PETG and ABS can be subjected to temperatures of about 150 degrees F before deforming. PLA will deform at about 140 degrees F.

Assembly

You may have received parts that require assembly. Many health care providers prefer to receive parts so they can more easily sanitize them before use.

Assembly is easy and documented in this short two minute video:

<https://youtu.be/qKqDUExM4Q8>

It is easy to take them apart. We have a short video on the process:

<https://youtu.be/F7FFPzGerME>

Support

If you have any issues or concerns please contact us at

faceshields@makernexus.org

More info

The latest information on the face shields, how to order more, and the design are available on our wiki: <http://makernexuswiki.com/index.php?title=Covid19>

Or our web site: <http://makernexus.org/covid19>